

Beer's law states that the amount of light absorbed by a chemical is *directly* related to the concentration of the chemical in a solution. A colorimeter (spectrophotometer) can provide useful qualitative and quantitative information about a chemical sample. This device measures the amount of light your samples absorb. The intensity of the light collected by a photocell is monitored as either an absorbance or a percent transmittance value. The value of this reading depends on the concentration of the coloured material in the solution.

### Question

How can Beer's law be used to determine an unknown concentration?

### Materials

computer system and interface  
colorimeter sensor  
cuvette  
5 test tubes (20 × 150 mm)  
test tube rack  
stirring rod  
transfer pipette  
100 mL of 0.40 M nickel(II) sulfate (NiSO<sub>4</sub>) solution  
100 mL distilled water

### Safety Precautions



- Nickel(II) sulfate is harmful if swallowed or inhaled. Nickel compounds are highly toxic by all routes of exposure and are irritants.
- Never pipette by mouth. Always use a pipette bulb or pipette pump.
- Be sure to dispose of materials properly.

### Procedure

1. Label the clean, dry test tubes #1 through #5.
2. Pipette 2, 4, 6, 8, and 10 mL of the NiSO<sub>4</sub> solution into test tubes #1 through #5 respectively.
3. Pipette 8, 6, 4, and 2 mL of the distilled water into test tubes #1 through #4 respectively.
4. Thoroughly mix the contents of test tubes 1 through 4 with the stirring rod, wiping it clean between stirrings. The following table outlines the volumes and concentrations that will be used:

Test tube #	Volume of NiSO <sub>4</sub> (mL)	Volume of water (mL)	Concentration (M)
1	2	8	0.08
2	4	6	0.16
3	6	4	0.24
4	8	2	0.32
5	10	0	0.40

**5.** Set up the computer system with the colorimeter sensor.

**6.** Set the recording of data to manual sampling control.

**7.** Calibrate the colorimeter sensor

**8.** Display the sensor with a table display of absorbance (A) and concentration (M).

**9.** Display the sensor with a digits display of absorbance.

**10.** Turn the wavelength setting of the colorimeter sensor to the RED (625 nm) LED colour.

**11.** Empty any material in the cuvette, and carefully rinse it with distilled water. Empty the water from the cuvette and dry it well.

**12.** Rinse the cuvette twice with approximately 1 mL amounts from test tube #1.

**13.** Fill the cuvette  $\frac{3}{4}$  full of solution from test tube #1 and carefully dry the outside of the cuvette.

**14.** Place the cuvette in the colorimeter and close the lid.

**15.** Once the absorbance digits display stabilizes, use the keyboard to save its value.

**16.** Continue testing each of the test tubes #2 through #5 by repeating steps 11 to 15.

**17.** Obtain the unknown concentration solution of nickel(II) sulfate from your teacher and repeat steps 11 to 15.

**18.** Record the absorbance for the unknown concentration solution.

**19.** Discard the solutions as directed by your teacher. Do not pour anything down the drain.

### **Analysis**

**1.** Create a graph display of absorbance (A) versus concentration (M).

**2.** Locate the absorbance value along the vertical axis for the unknown solution, and trace a horizontal line to the plotted line.

**3.** What is the corresponding horizontal value directly below the plotted line?

### **Conclusions**

**4.** Write a conclusion to explain how your experimental observations supported your theoretical calculations.

### **Applications**


**5.** How can the use of Beer's law and spectroscopy be useful in industries?

**6.** What would be the limitations of using this technique in industries?

## Teacher Information

### Sample Data – Answer Sheet

#### Observations



**Graph 1**

* Light Absorbance Run #1	
Absorbance (A)	Concentration (M)
17.440	0.080
93.221	0.160
141.358	0.240
188.216	0.320
228.207	0.400

**Table 1**

#### Answers

**5.** Concentrations of unknown solutions can be measured quickly and easily once the Beer's law graph is determined.

**6.** It is limited to solutions that are not completely colourless or opaque.